



# READY LABEL

## Label Design Simplified



### WHY READY LABEL?

Design and print custom labels using our simple but flexible, touch-optimized label design software. Use Ready Label in a stand-alone print station or integrate it into a fully automated system. However you plan to print your labels, Ready Label will work for you.

### PURCHASING OPTIONS:

<b>PRINT STATION:</b>	Includes tablet with software installed, docking station, bluetooth keyboard, mouse, CAB label printer and (1) each rolls of labels and ribbon (standard size).
<b>PLUG AND PRINT:</b>	Tablet with software installed, docking station, bluetooth keyboard and mouse
<b>SOFTWARE:</b>	Individual software license



### FEATURES SET:

- Import data from Excel and .CSV files
- Wide variety of 1D and 2D barcodes: code 128, code 39, Datamatrix, QR, etc.
- Generate custom alpha-numeric sequences, and time and date options
- User selectable fonts
- Import custom images such as company or product logos
- Fully customizable print layout
- Can be integrated with Green Button Go™ and other automation scheduling software

### COMPATIBLE PRINTERS:

- CAB Printers
  - Zebra Printers
  - Brother Printers
  - Cedrex Biolabel
  - Kbiosystems Gecko
  - Agilent VCode
- Drivers are available for other third party automated labelers*

### RECOMMENDED PC REQUIRMENTS:

- Intel Core™ i3 or Intel Atom equivalent processor
- 2 GB RAM
- Compatible with Windows XP, Vista, 7 and 8
- For any further questions regarding system requirements and capabilities contact your local Biosero representative